


# PROVOX LaryButton



1


2


3


3a, a-1


a-2


3b, b 1


b-2


b-3


b-4


4a

4b


4c


5a


5b


5c


5d

5e

5f

5g

## **Prescription information**

**CAUTION:** United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

## **Disclaimer**

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

## Contents

EN - ENGLISH .....	9
DE - DEUTSCH .....	15
NL - NEDERLANDS .....	23
FR - FRANÇAIS .....	30
IT - ITALIANO .....	37
ES - ESPAÑOL .....	45
PT - PORTUGUÊS .....	52
SV - SVENSKA .....	59
DA - DANSK .....	66
NO - NORSK .....	73
FI - SUOMI .....	79
IS - ÍSLENSKA .....	87
ET - EESTI .....	94
LT - LIETUVIŲ K. ....	100
CS - ČESKY .....	107
HU - MAGYAR .....	114
SK - SLOVENČINA .....	121
SL - SLOVENŠČINA .....	128
PL - POLSKI .....	135
RO - ROMÂNĂ .....	142
HR - HRVATSKI .....	149
SR - SRPSKI .....	156
EL - ΕΛΛΗΝΙΚΑ .....	163
TR - TÜRKÇE .....	171
HY - ՀԱՅԵՐԵՆ .....	178
RU - РУССКИЙ .....	186
MS - BAHASA MELAYU .....	194
JA - 日本語 .....	201
KO - 한국어 .....	211
ZH-TW - 繁體中文 .....	219
HE - עברית .....	233
AR - العربية .....	238

## 1. Descriptive information

### 1.1 Intended use

The Provox LaryButton is a self-retaining holder for devices in the Provox HME System intended for vocal and pulmonary rehabilitation after total laryngectomy.

For patients with a shrinking tracheostoma it is also used to maintain the tracheostoma for breathing.

The Provox LaryButton is intended for single patient use.

### 1.2 CONTRAINDICATIONS

Provox LaryButton is not intended to be used by patients that:

- are under any form of mechanical ventilation.
- have damaged tracheostoma tissue.

### 1.3 Description of the device

The Provox LaryButton is a self-retaining tracheostoma button made of medical grade silicone rubber. The purpose of the device is to create a comfortable and airtight fit between the Provox LaryButton and the tracheostoma, and also to provide attachment for devices from the Provox HME System. The device is delivered single packed, non-sterile and ready to use.

The different parts of the Provox LaryButton are (Fig. 1):

- a) Retention Collar
- b) Shaft
- c) Shield (conical)
- d) HME and Accessory Holder
- e) Wings

### 1.4 WARNINGS

- Always use a Provox LaryButton of an appropriate size. If the Provox LaryButton is too wide or too long it may cause tissue damage or irritation. If it is too long it may also cause dislodgement of the voice prosthesis or XtraFlange. Furthermore, insertion or removal of the HME Cassette as well as pressing on the HME Cassette to speak may exert pressure on the voice prosthesis if the Provox LaryButton is too long.
- Always make sure that the HME Holder part is large enough in relation to the size of the patient's stoma. If the product is pushed too far into trachea it may obstruct breathing.

- **DO NOT** cut off the wings or otherwise mechanically modify the Provox LaryButton since it will change the outer diameter and/or the mechanical stability of the product. This may result in the LaryButton being pushed into the trachea, which may obstruct breathing. A modified LaryButton may also cause irritation and/or bleeding.
- **DO NOT** fenestrate the LaryButton as the fenestration might cause dislodgement of the voice prosthesis or XtraFlange.
- **DO NOT** reuse the Provox LaryButton between patients. The device is intended for single patient use only. Reuse between patients may cause cross contamination.
- Always insert and remove the Provox LaryButton according to the instructions below.  
**DO NOT** use force and do not fold the part of the LaryButton that holds the HME Cassette during insertion. If the product is pushed too far into trachea it may obstruct breathing and the use of force during insertion may cause bleeding.
- If used, **DO NOT** push the tip of the hemostats or forceps through the Provox LaryButton (Fig. b-2). If the instrument tip is protruding from the button during insertion, the tissue of the trachea could be harmed and/or the voice prosthesis could become dislodged and subsequently be aspirated or ingested.
- **DO NOT** use lubricating gels if the patient has allergies related to these substances.
- Always remove the device from the stoma before cleaning. **DO NOT** at any time attempt to clean the device while inside the stoma since this could cause tissue damage and dislodgement of voice prosthesis.
- Make sure that the patient has been trained in the use of the device, and has demonstrated the ability to understand and consistently follow the Instructions for Use without clinician supervision.

## 1.5 PRECAUTIONS

- **DO NOT** use a dirty or contaminated Provox LaryButton. Clean and disinfect the device according to the cleaning and disinfection instructions below.
- **DO NOT** use toothed or sharp instruments at any time during insertion since they could damage the LaryButton and cause tissue damage.
- **DO NOT** use a damaged device since this might cause aspiration of small parts or cause formation of granulation tissue.
- Only use water-soluble lubricants. Oil-based lubricants (e.g., Vaseline) may weaken, damage, or destroy the product and should not be used at any time.
- Carefully monitor the tracheostoma tissue during radiation therapy and stop using the device if the stoma gets irritated or starts to bleed.

- Patients with bleeding disorders or patients on anticoagulants should not use the device if it causes bleeding that re-occurs.
- Stop using the device if granulation tissue forms in the tracheostoma.

## 2. Instructions for use

### 2.1 Preparation

#### Choosing correct size

**Diameter:** To achieve an airtight fit and good retention the diameter of the Retention Collar of the LaryButton should exceed the diameter of the tracheostoma entrance. This is especially important for hands-free speech. Measure the size of the tracheostoma (largest diameter) with a ruler and fit the Provox LaryButton accordingly. In irregularly shaped tracheostomas, it may be difficult to achieve an airtight fit. Selection of a too small diameter will lead to a too loose fit; the device may fall out or air may escape when speaking. Selection of a too wide LaryButton may irritate or damage the tracheostoma.

**Length:** In most cases the shorter size where the Retention Collar of the Provox LaryButton fits right behind the edge of the tracheostoma will provide a correct fit. However if the device is too short it may fall out of the tracheostoma.

To help to choose what size(s) of LaryButton to prescribe to the patient, the clinician can use a Provox LaryButton Sizer Kit. The Sizer Kit contains samples (“sizers”) of commercially available Provox LaryButtons.

**CAUTION:** During fitting always make sure that the LaryButton does not reach the tracheal flange of the voice prosthesis.

Optimal fit and performance of the Provox LaryButton is achieved if the patient has a so-called “tracheostoma lip” that is situated uniformly around the tracheostoma.

#### Preparing for insertion

Before first use check that you have the correct size and that the package is not damaged or opened. Do not use otherwise.

Make sure that your hands are thoroughly clean and dry before placing them in the area of the stoma or handling the Provox LaryButton and/or accessories.

Have a mirror and sufficient light to illuminate the stoma.

Always carefully inspect the Provox LaryButton BEFORE EACH USE (i.e. before insertion), and ensure that it is undamaged and shows no signs of tears, cracks, or crusts. If you detect any damage, do not use the product, and obtain a replacement.

Before insertion, rinse the device with drinking water to remove any disinfectant residues (Fig. 2). Make sure that you only insert clean devices.

In a hospital setting it is better to use sterile water rather than drinking water.

If needed, lightly lubricate the Retention Collar of the device with some water-soluble lubricant. The following water-soluble lubricants are recommended for use with the Provox LaryButton:

- KY Jelly<sup>®</sup>
- SurgiLube<sup>®</sup>
- Xylocaine<sup>®</sup>/(lidocaine)

**CAUTION:** Only use water-soluble lubricants. Oil-based lubricants (e.g., Vaseline) may weaken damage or destroy the product and should not be used at any time.

Make sure that you do not lubricate the HME Holder, HME Cassette or any accessory which is to be held by the Provox LaryButton, as this may lead to accidental detachment.

## 2.2 Operating instruction

### Insertion (Fig. 3):

#### Alternative 3a:

1. Fold the Retention Collar with your fingers (Fig. 3) and hold it firmly folded.
2. Insert the Retention Collar into the stoma and release it (Fig. a-1 and a-2).

#### Alternative 3b:

1. Fold the Retention Collar with your fingers (Fig. 3).
2. Use a pair of blunt forceps or a lockable blunt hemostat from the inside to keep it folded (Fig. b-1 and b-2).

**CAUTION: DO NOT** push the tip of the hemostats or forceps through the Provox LaryButton (Fig. b-2). If the instrument tip is protruding from the button during insertion, the tissue of the trachea could be harmed and/or the voice prosthesis could become dislodged and subsequently be aspirated or ingested.

3. Insert the Retention Collar into the stoma (Fig. b-3).
4. Release the forceps to unfold the button (Fig. b-4).

Check the correct fit by slightly pulling and/or rotating the button.

Gently insert the Provox HME System component into the HME Holder.

**NOTE:** Some patients might experience more coughing after inserting the LaryButton and when touching the device to speak. This usually diminishes when the user gets used to the device. If discomfort occurs the LaryButton should be removed.

### Device removal

1. Hold the Provox LaryButton in place and gently remove the Provox HME

System component from the HME Holder (Fig. 4a).

2. Grasp the Provox LaryButton and gently turn it (about half a turn) to allow the Retention Collar to collapse thereby reducing its width (Fig. 4b).
3. Remove the device by gently tilting it out of the stoma (Fig. 4c). Alternatively, the device can be removed by pushing with a finger on the shaft from one direction and tilt the device out.

**NOTE:** The insertion or removal procedures may occasionally cause some slight bleeding, irritation or coughing. Inform your clinician if the bleeding persists.

### 2.3 Cleaning and disinfection

**CAUTION:** Always remove the device from the stoma before cleaning.

**DO NOT** at any time attempt to clean the device while inside the stoma since this could cause injury.

Clean the Provox LaryButton after each use, or at least daily with the Provox TubeBrush (Fig. 5a-5g).

Do not use any water other than drinking water to clean and rinse the device.

**CAUTION:** Be careful not to burn your finger in the warm cleaning water.

Disinfect the Provox LaryButton at least once a day (Fig. 5e). with one of the following methods:

- Ethanol 70% for 10 minutes
- Isopropylalcohol 70% for 10 minutes
- Hydrogenperoxide 3% for 60 minutes

**CAUTION: DO NOT** use the device until it is completely dry (Fig. 5f). Inhalation of disinfectant fumes can cause severe coughing and irritation of the airways.

If the Provox LaryButton looks dirty or has air dried in an area with a risk of contamination the device should be both cleaned and disinfected before use. A risk of contamination could exist if the device has been dropped on the floor or if it has been in contact with a pet, someone with a respiratory infection, or any other gross contamination.

During hospitalization it is important to both clean and disinfect the Provox LaryButton after use but also before use since there is an increased risk of infection and contamination. In a hospital it is better to use sterile water for cleaning and rinsing, rather than drinking water.

**CAUTION: DO NOT** clean or disinfect by any other method than described above as this might cause product damage and patient injury.

### 2.4 Storage instructions

When not in use, clean and disinfect the device as described above and then

store the Provox LaryButton in a clean and dry container at room temperature. Protect from direct sunlight.

## 2.5 Device lifetime

The Provox LaryButton may be used for a maximum period of 6 months. If the device shows any signs of damage it shall be replaced earlier.

## 2.6 Disposal

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

## 2.7 Accessories

**CAUTION:** Use only original components. Other accessories may cause product damage and/or malfunction.

Some of the devices that can be kept in place by the Provox LaryButton are:

- Provox HME Cassettes (Heat and Moisture Exchangers)
- Provox FreeHands HME, for speaking without manual stoma occlusion.
- Provox ShowerAid that can help to prevent water from entering the LaryButton during showering.

If needed a Provox TubeHolder or Provox LaryClips can be used for additional support to prevent the Provox LaryButton from falling out of the stoma.

## 3. Troubleshooting information

### Bleeding of the stoma

Stop using the device. After the tracheostoma is healed try gradually increasing the duration of using the device. If even short periods of use cause bleeding or if the bleeding persists after removing the device, stop using it and consult your clinician.

### Rapid shrinking of the stoma

If you are wearing the Provox LaryButton because of a rapidly shrinking tracheostoma, discuss with your clinician the possibility of having two devices (of the same size). This would allow you to remove the Provox LaryButton from your stoma and replace it immediately with a clean one. Cleaning of the removed device can then be done afterwards.

### Stoma enlarges during use

If your tracheostoma enlarges after using the device for some hours, ask your clinician if it is possible for you to use two devices (one smaller and one larger diameter). This allows you to use the larger size when the tracheostoma has enlarged, and the smaller size to allow the tracheostoma to shrink. It may be necessary to discontinue use of the button during the night, so that the

tracheostoma can shrink. In that case, discuss alternative options to attach the HME Cassette with your clinician so that you can continue HME Cassette use during the night.

## **4. Additional information**

### **4.1 Compatibility with MRI Examination**

**MR-Safe:** This device does not contain any metallic elements and has no potential for interaction with the MRI field.

### **4.2 User assistance information**

For additional help or information please see back cover of this Instructions for Use for contact information.

## **5. Reporting**

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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Atos and the Atos logo are trademarks of Coloplast A/S.
© 2025 Coloplast A/S. All rights reserved.
**Atos Medical AB**