

## Tracoe Twist Plus – REF 523

Double fenestrated, with 15 mm connector  
(blue locking ring)

### Product Description

The Tracoe Twist Plus tracheostomy tubes are used as ventilation and therapy tubes with inner cannulas in hospitals, rehabilitation, and community settings. A special feature is the curved neck flange that moves in both the vertical and horizontal axes, following the movements of the patient's head. The thin tube walls of inner and outer cannulas allow for a good inside - outside diameter ratio. Tracoe Twist and Twist Plus tracheostomy tubes are made of polyurethane (PU). The low-pressure cuff is made of medical grade PVC.



### Scope of Supply

- 3 inner cannulas, double fenestrated with 15 mm connector

### General Information

Product line	Twist Plus
Legal manufacturer	Tracoe Medical GmbH Reichelsheimer Str. 1/3 55268 Nieder-Olm, Germany
Country of origin	Germany
Colour Coding	n. A.
Handling	
Shelf life	5 years, from date of production
Maximum Useful Life	5 weeks, from first use

### Packaging

Packaging unit	3 units
Weight	See Table
Dimensions (W x H x L)	93 x 200 x 60
Material outer packaging:	Cardboard box
Material sterile packaging:	Sterile bag with medical paper and transparent heatsealable PET/PP foil

### Material Specifications

Material Outer cannula	Polyurethane
Material Neck Flange	n. A.
Material Cuff	n. A.
Biocompatibility	EN ISO 10993-1



### Storage



### Certification

Product classification acc. MDR	IIb
CE Marking	
GMDN	35404
Quality system certified to	EN ISO 13485
FDA registration	Clearance Number K961449

## Product Range

Size	Order Number	ID-IC	BA	Gross Weight	Global Trade Item Number
		mm	°	g	
07	REF 523-07	7,0	100	42	4035324018277
08	REF 523-08	8,0	100	41	4035324018284
09	REF 523-09	9,0	100	42	4035324018291
10	REF 523-10	10,0	100	43	4035324018307

## Legend

- ID-IC** Inside diameter at bottom of inner cannula
- BA** Bending angle



Tracoe Medical GmbH  
Reichelsheimer Str. 1 / 3  
55268 Nieder-Olm  
Germany



Atos Medical AB  
Hyllie Boulevard 17  
215 32 Malmö  
Sweden  
[www.atosmedical.com](http://www.atosmedical.com)