

Tracoe Technical Equipment – REF 721 Tracoe Cuff Pressure Monitor Sensitive

Product Description



The handheld cuff pressure monitor is used to fill and monitor the pressure of High-Volume Low-Pressure cuffs of tracheostomy and endotracheal tubes. It has an Easy-to-read scale from 0 to 120 cmH₂O.





The usual pressure range between 20 and 30 cmH₂O is marked in green and allows a precise, reproducible pressure indication.

Scope of Supply


- Tracoe Cuff Pressure Monitor Sensitive 1pc
- 100 cm Connection Tube 1pc



General Information	
Product line	Technical Equipment
Legal manufacturer	VBM Medizintechnik GmbH Einsteinstraße 1 72172 Sulz a.N. Germany
Country of origin	Germany
Handling	 
Shelf life	8 years for CPM 5 years for Connection Tube

Material Specifications	
 	
Product specifications	
Net weight CPM	156 g
Storage	
 	

Packaging	
Packaging unit	1 piece
Dimensions (W x H x L)	210 x 70 x 81 mm
Gross weight	234 g
Material outer packaging:	Cardboardbox

Certification	
Product classification acc. MDR	I
CE marking	
Quality system certified to	EN ISO 13485
FDA registration	Reg. No. 8685750



Tracoe Medical GmbH,
Reichelsheimer Str. 1 / 3
55268 Nieder-Olm
Germany



Breathing-Speaking-Living

Atos Medical AB
Hyllie Boulevard 17
215 32 Malmö
Sweden
www.atosmedical.com