

Tracoe Technical Equipment – REF 730-5 Tracoe Smart Cuff Manager

Product Description

The category Technical Equipment refers to medical devices that allow the pressure in the cuffs of tracheostomy and endotracheal tubes to be set or independently regulated with great precision.

The Tracoe Smart Cuff Manager monitors and regulates the internal pressure of High-Volume Low-Pressure Cuffs) of tracheostomy tubes as well as endotracheal tubes.



Scope of Supply

- Tracoe Smart Cuff Manager
- Tracoe connection tube – REF 702

General Information

Product line	Technical Equipment
Legal manufacturer	Tracoe Medical GmbH Reichelsheimer Straße. 1/3. 55268 Nieder-Olm Germany
Country of origin	Germany
Handling	
Shelf life	2 years
Maximum Useful Life	29 days

Packaging

Packaging unit	5 pieces
Dimensions (W x H x L)	43 x 80 x 161 mm
Weight	312 g
Material outer packaging:	Cardboardbox
Material primary packaging:	Sterile bag with medical paper and transparent heatsealable PET/PP foil

Material Specifications

Material Housing	Acrylonitrile buta-diene styrene (ABS)
Material Ballon	Termoplastic elastomer (TPE)
Biocompatibility	ISO 10993-5 and 10993-12



Product specifications

Weight	16 g
Dimensions (L x H)	55 x 77 mm

Storage



Certification

Product classification acc. MDR	IIa
CE Marking	0123
GMDN	35401
Quality system certified to	EN ISO 13485
FDA registration	Clearance Number K152865